


Product Name: Recombinant Human PALMD Protein
Catalog #: PEH8888



Summary

Name	PALMD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NP74
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ANIPDLPSY IPSRLRKEIN EEKEDDEQNR KALYAMEIKV EKDLKTGEST VLSSIPLPSD DFKGTGIKVV DDGQKSVYAV SSNHSAAAYNG TDGLAPVEVE ELLRQASERN SKSPTEYHEP VYANPFYRPT TPQRETVTPG PNFQERIKIK TNGLGIGVNE SIHNMGNGLS EERGNNFNHI SPIPPVPHPR SVIQQAEELK HTPQKRLMTP WEESNVMQDK DAPSPKPRLS PRETIFGKSE HQNSSPTCQE

image



Background

Alternative Names Palmdelphin (Paralemmin-like protein)

Note

Product Name: Recombinant Human PALMD Protein
Catalog #: PEH8888



For Research Use Only , Not for Diagnostic Use