


Product Name: Recombinant Human PALM2 Protein
Catalog #: PEH8889



Summary

Name	PALM2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IXS6
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ERLQAIAEKR KRQTEIEGKR QQLDEQILLL QHSKSKVLRE KWLLQGIPAG TAESEEARRR QSEDEFVRVK QLEDNIQRLE QEIQTLESEE SQISAKEQII LEKLEKETEKS FKDFQKGFSS TDGAVYAMEI NVEKDKQTGE TKILSTSTIG PEGVHQKGVK VYDDGTVVY EVRSGGTVE NGVHKLSTKD VEELIQKAGQ SSLGGGHVSE RTVIADGSL S HPKEHMLCKE AKLEMVHKS R KDHSSGNPGQ QAQAPSAAGP EANLDQPVTM IFMGYQNIED EETKVKLVG DETIKAELV L IDEDEKSLR EKTVDVSTI DGNAAEVSG R PVSDTTEPS

image



Background

Alternative Names Paralemmin-2

Product Name: Recombinant Human PALM2 Protein
Catalog #: PEH8889



Note

For Research Use Only , Not for Diagnostic Use