


Summary

Name	RIMB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95153
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GDPGAMEPWA LPTWHSWTPG RGGEPSSAAP SIADTPPAAL QLQELRSEES SKPKGDGSSR PVGGTDPEGA EACLPSLGQQ ASSSGPACQR PEDEEVEAFL KAKLNMSFGD RPNLELLRAL GELRQRCAIL KEENQMLRKS SFPETEEKVR RLKRKNAELA VIAKRLEERA RKLQETNLRV VSAPLPRPGT SLELCRKALA RQRARDLSET ASALLAKDKQ IAALQRECRE LQARLTLVGK EGPQWLHVRD FDRLRESQR EVLRLQRQIA LRNQRETLPL PPSWPPGPAL QARAGAPAPG APGEATPQED ADNLPVILGE PEKEQRVQQL ESELSKKRKK CESLEQEARK KQRRCEELEL QLRQAQNENA RLVEENSRLS GRATEKEQVE

image



Background

Alternative Names Peripheral-type benzodiazepine receptor-associated protein 1 (PRAX-1)

Product Name: Recombinant Human TSPOAP1 Protein
Catalog #: PEH8891



(Peripheral benzodiazepine receptor-interacting protein) (PBR-IP) (RIMS-binding protein 1) (RIM-BP1) (TSPO-associated protein 1)

Note

For Research Use Only , Not for Diagnostic Use