


**Product Name: Recombinant Human NUDT7 Protein**  
**Catalog #: PEH8892**



## Summary

<b>Name</b>	NUDT7_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	POC024
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	28 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	MSRLGLPEEP VRNLLDDAK ARLRKYDIGG KYSHLPYNKY SVLLPLVAKE GKLHLLFTVR SEKLRRAPGE VCFPGGKRDPTDMDDAATAL REAQEEVGLR PHQVEVVCCLVPLIDDTL ITPFVGLIDH NFQAQNPAAE VKDVFLVPLA YFLHPQVHDQ HYVTRLGHRF INHIFEYTNP EDGVTYQIKG MTANLAVLVA FIILEKKPTF EVQFNLNDVL ASSEELFLKV HKKATSRL

## image



## Background

**Alternative Names** Peroxisomal coenzyme A diphosphatase NUDT7 (EC 3.6.1.-) (Nucleoside diphosphate-linked moiety X motif 7) (Nudix motif 7)

## Note

**Product Name: Recombinant Human NUDT7 Protein**  
**Catalog #: PEH8892**



---

For Research Use Only , Not for Diagnostic Use