


Product Name: Recombinant Human PHF24 Protein
Catalog #: PEH8894



Summary

Name	PHF24_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UPV7
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MGVLM SKRQT VEQVQKVS LA VSAFKDGLRD RPSIRRTGEL PGSRRGTVEG SVQEVQEEKE AEAGTSVVQE ESSAGRAAWE RLRDGRGVEP EEFDRTSRFT PPAFIRPTRK LDDDKPPEIC LEPREPVVND EMCDVCEVWT AESLFPCRVC TRVFHDGCLR RMGYIQGDSA AEVTEMAHTE TGWSCHYCDN INLLLTEEM YSLTETFQRC KVIPDCSLTL EDFLRYRHQA AKRGDRDRAL SEEQEEQAAR QFAALDPEHR GHIEWPDFLS HESLLLLQQL RPQNSLLRLL TVKERERARA AFLARGSGST VSEAECRRAQ HSWFCKRFPE APSCSVSISH VGPIADSSPA SSSSKSQDKT LLPTEQESRF VDWPTFLQEN VLYILAARPN SAAIHLKPPG

image



Background

Alternative Names PHD finger protein 24

Product Name: Recombinant Human PHF24 Protein
Catalog #: PEH8894



Note

For Research Use Only , Not for Diagnostic Use