


Product Name: Recombinant Human PCARE Protein
Catalog #: PEH8897



Summary

Name	PCARE_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A6NGG8
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IPLLVKNSTC YDAGEGLAEE QPSRRNQTT AKGLCQLMGD PASGKRKDME GLIPGKTSS SQLNKSQSHM AKDIPFKTQG SHGSQGADFS GDESEESSTQ DTSKWKRKTAK CHTSSTQSHC YQTIHPAHEP EGKVDPEPL VKAHQQAYTY

image



Background

Alternative Names Photoreceptor cilium actin regulator

Note

For Research Use Only , Not for Diagnostic Use