


Product Name: Recombinant Human PLEKHA3 Protein
Catalog #: PEH8898



Summary

Name	PKHA3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HB20
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PDPLVSPVSP SPVQMMKRSV SHPGSCSSER SSHSIKEPVS TLHRLSQRRR RTYSDTDCS DIPLEDPDRP VHCKNTLNG DLASATIP EE SRLMAKKQSE

image



Background

Alternative Names	Pleckstrin homology domain-containing family A member 3 (PH domain-containing family A member 3) (Phosphatidylinositol-four-phosphate adapter protein 1) (FAPP-1) (Phosphoinositol 4-phosphate adapter protein 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use