


Product Name: Recombinant Human DDX10 Protein
Catalog #: PEH8902



Summary

Name	DDX10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13206
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TKELVRSQAD KVI EPRAPSL TNDEVEEFRA YFNEKMSILQ KGGKRLEGTE HRQDNDTGNE EQEEEEEDDEE EMEEKLAKAK GSQAPSLPNT SEAQKIKEVP TQFLDRDEEE EDADFLKVKR HNVFGLDLKD EKT LQKKEPS KSSIKKKMTK VAEAKKVMKR NFKVNNKITF TDEGELVQQW PQMQKSAIKD

image



Background

Alternative Names Probable ATP-dependent RNA helicase DDX10 (EC 3.6.4.13) (DEAD box protein 10)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human DDX10 Protein
Catalog #: PEH8902

