

Product Name: Recombinant Human MKRN2 Protein
Catalog #: PEH8906



Summary

Name	MKRN2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H000
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ENPIIKSCPE CRVISEFVIP SVYWVEDQNK KNELIEAFKQ GMGKKACKYF EQGKGTCPPG SKCLYRHAYP DGRLAPEKP RKQLSSQGTV RFFNSVRLWD FIENRESRHV PNNEDVDMTE LGDLFMHLSG

image



Background

Alternative Names	Probable E3 ubiquitin-protein ligase makorin-2 (EC 2.3.2.27) (RING finger protein 62) (RING-type E3 ubiquitin transferase makorin-2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use