



Product Name: Recombinant Human LARS2 Protein
Catalog #: PEH8909



Summary

Name	SYLM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15031
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AMQDALADLP EWYGIKGMQA HWIGDCVGCH LDFTLKVHGQ ATGEKLTAYT ATPEAIYGTS HVAISPSHRL LHGHSSLKEA LRMALVPGKD CLTPVMAVNM LTQQEVPVVI LAKADLEGLS DSKIGIPSTS SEDTILAQTL GLAYSEVIET LPDGTERLSS SAEFTGMTRQ DAFLALTQKA RGKRVGGDVT

image



Background

Alternative Names	Probable leucine--tRNA ligase, mitochondrial (EC 6.1.1.4) (Leucyl-tRNA synthetase) (LeuRS)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human LARS2 Protein
Catalog #: PEH8909

