


Product Name: Recombinant Human HYCC2 Protein
Catalog #: PEH8911



Summary

Name	F126B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IXS8
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SVNLNDADEG FSSGASLSSQ PIGTKPSSSS QRGSLRKVAT GRSKDKETA SAIKSSESPR DSVVRKQYVQ QPTDLSVDSV ELTPMKKHL LPAQQVVPKI NSLSLIRTAS ASSSKSFDYV NGSQASTSIG VGTEGGTNLA ANNANRYSTV SLQEDRLGQA GEGKELLSPG APLTKQSRSP

image



Background

Alternative Names Protein FAM126B

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human HYCC2 Protein
Catalog #: PEH8911

