


Product Name: Recombinant Human FAM83A Protein
Catalog #: PEH8913



Summary

Name	FA83A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86UY5
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DFLSSVEAQY IQAQAREPPC PPDTLGGAEA GPKGLDSSSL QSGTYFPVAS EGSEPALLHS WASAEKPYLK EKSSATVYFQ TVKHNNIRDL VRRICITRTSQ VLVILMDVFT DVEIFCDILE AANKRGVFVC VLLDQGGVKL FQEMCDKVQI SDSHLKNISI RSVEGEIYCA KSGRKFAGQI REKFIISDWR FVLSGYSYFT WLCGHVHRNI LSKFTGQAVE LFDEEFRHLY ASSKPVMGLK SPRLVAPVPP

image



Background

Alternative Names Protein FAM83A (Tumor antigen BJ-TSA-9) (Tumor-specific gene expressed in prostate protein)

Note

Product Name: Recombinant Human FAM83A Protein
Catalog #: PEH8913



For Research Use Only , Not for Diagnostic Use