


Product Name: Recombinant Human FEM1C Protein
Catalog #: PEH8915



Summary

Name	FEM1C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96JP0
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ARDGKLRLLT KLLASKSKEE VSSLISEKTN GATPLLMAAR YGHLDMVEFL LEQCSASIEV GGSVNFDFGET IEGAPPLWAA SAAGHLKVVQ SLLNHGASVN NTTLTNSTPL RAACFDGHLE IVKYLVEHKA DLEVSNRHGH TCLMISCYKG HKEIAQYLLE KGADVNRKSV KGNTALHDCA ESGSLDIMKM LLMYCAKMEK DGYGMTPLLS ASVTGHTNIV DFLTHHAQTS KTERINALEL LGATFVDKRR DLLGALKYWK KAMNMRYSDR TNIISKVPVQ TLIMAYDYAK EVNSAAEELEG LIADPDEMRRM

image



Background

Alternative Names Protein fem-1 homolog C (FEM1c) (FEM1-gamma)

Product Name: Recombinant Human FEM1C Protein
Catalog #: PEH8915



Note

For Research Use Only , Not for Diagnostic Use