



Product Name: Recombinant Human FRY Protein
Catalog #: PEH8917



Summary

Name	FRY_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5TBA9
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AVTRSTSSTS SGSNSNVLP VSWKRPQYSQ KRTKEKLVHV LSLCGQEVGL SKNPSVIFSS CGDLLEHQ TSLVSEEDGA REQENMDDTN SEQQFRVFRD FDFLDVELED GEGESMDNFN WGVRRRSLDS LDKCDMQILE ERQLSGSTPS LNKMHEDSD ESSEEDLTA

image



Background

Alternative Names Protein furry homolog

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FRY Protein
Catalog #: PEH8917

