



Product Name: Recombinant Human KATNIP Protein
Catalog #: PEH8919



Summary

Name	KATIP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O60303
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	APADHSILVD QKNEKSEQL EAMNAHSEES KGTHEMAGAS GDKELGLGCS PPAETLADAK LSSQGNVSGK RKNSTNCRKD SLSQLEEYLR LSAVPTSMGD MPSAPATSPP VKCPPVHEEP SLIQQLENLM GRKICEPPGK TPSWLQPSPT GKDRKQGGRK PKPLWLSPEK PLAWKGR LPS DDVIGEGPGE TEARDKGLRH EPGWGTSRSV NTKERPQRAT TKVHSDDSDI FNQPPNRERP ASGRRGSRKD AGSSSHGDDQ PASREDTWSS RTPSRSRWRS EQEHTLHESW SLSAFDRSH RGRISNTELP GDILDELLQQ KSSRHSDLPP SKKGEPGLS RGQDGYSGET

image



Background

Alternative Names Katanin-interacting protein

Product Name: Recombinant Human KATNIP Protein
Catalog #: PEH8919



Note

For Research Use Only , Not for Diagnostic Use