



Product Name: Recombinant Human KRI1 Protein
Catalog #: PEH8920



Summary

Name	KRI1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N9T8
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RLKDRYGDRD SSSDSSSESD SSERVEFDP QQERDFYKTL SLLKKKDPRI YQKDATFYNR TASSDSEED PEALEKQKKV RPMYLKDYER KVILEKAGKY VDEENS DGET SNHRLQETSS QSYVEEQQL KESFRAFVED SEDEDGAGEG GSSLLQKRAK TRQEKAQEEA DYIEWLKGQK EIRNPDSLKE LTHLKEYWND PELDEGERFL

image



Background

Alternative Names Protein KRI1 homolog

Note

Product Name: Recombinant Human KRI1 Protein
Catalog #: PEH8920



For Research Use Only , Not for Diagnostic Use