

Product Name: Recombinant Human PRR16 Protein
Catalog #: PEH8921



Summary

Name	LARGN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q569H4
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSCPAEGPPA ASKTKVKEQI KIIVEDLELV LGDLKDVAKE LKEVVDQIDT LTSDLQLEDE MTDSSKDTL NSSSSGTTAS SLEKIKVQAN APLIKPPAHP SAILTVLRKP NPPPPPPRLT PVKCEDPKRV VPTANPVKTN GTLLRNGGLP GGPNKIPNGD ICCIPNSNLD KAPVQLLMHR PEKDRCPQAG PRERVRFNEK VQYHGYCPDC DTRYNIKNRE VHLHSEPVHP

image



Background

Alternative Names Protein Largin (Mesenchymal stem cell protein DSC54) (Proline-rich protein 16)

Note

Product Name: Recombinant Human PRR16 Protein
Catalog #: PEH8921



For Research Use Only , Not for Diagnostic Use