


Product Name: Recombinant Human C17orf97 Protein
Catalog #: PEH8922



Summary

Name	LIAT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZQX7
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RAESRRLVGI GPRAPPGRVG LQPSGRLLDRR GGAGTMGYKD NDGEEEEEREG GAAGPRGSRL PPITGGASEL AKRKVKKKKR KKKTKGSGKG DDKHQSQLK SQPLSSSFHD ILSPCKERGP KPEHRQSKVE KKHLPDSST

image



Background

Alternative Names Protein LIAT1 (Ligand of ATE1 protein)

Note

For Research Use Only , Not for Diagnostic Use