


Product Name: Recombinant Human C11orf65 Protein
Catalog #: PEH8924



Summary

Name	MFI_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NCR3
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KQDKAARVIQ QAWKSFLNVA IFQHFSLID LRRQGEPRQI VKYINPKEAE LLDAAAGIHV RFRLGGVKFP PDIYYKIFTH RPIEDLCANS PRNYAKLPAK HTSHNKNDHL QEEDHSGWYH RIENNGWRPV SDTFWLSTDG MVVEDKKESE FHFSKLKRRQ DLEKKRKLK IEWMRQMYYS GSLEAKSTHH ETLGLIHTAT KGLIRAFEDG GIDSVMWEV

image



Background

Alternative Names Protein MFI (Mitochondrial fission factor interactor)

Note

Product Name: Recombinant Human C11orf65 Protein
Catalog #: PEH8924



For Research Use Only , Not for Diagnostic Use