


**Product Name: Recombinant Human C17orf75 Protein**  
**Catalog #: PEH8928**



## Summary

<b>Name</b>	NJMU_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9HAS0
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	39 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	GDEKELESSE EGGSAEERRL EPPSSSHYCL YSYRGSRLAQ QRGDSEDGSP SGTNAETPSG DDFSLSLADT NLPSEVEPEL RSFIKRLSR GAVFEG LGNV ASVELKIPGY RVGCYCLFQ NEKLLPETVT IDSERNPSEY VVCFLGGSEK GLELFRLELD KYIQGLKNNM NCEARGLESH IKSYLSSWFE DVVCPPIQRVV LLFQEKLTF LHAALSYPV EVKESDEKTK RDINRFLSVA SLQGLIHEGT MTSLCMAMTE EQHKSVVIDC SSSQPQFCNA GSNRFCDWM QAFLNGAKGG NPFLFRQVLE NFKLKAIQDT NNLKRFIRQA

## image



## Background

**Alternative Names** Protein Njmu-R1

**Product Name: Recombinant Human C17orf75 Protein**  
**Catalog #: PEH8928**



---

**Note**

For Research Use Only , Not for Diagnostic Use