


Product Name: Recombinant Human PPP1R15B Protein
Catalog #: PEH8930



Summary

Name	PR15B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q5SWA1
Host	E.coli
Species	Human
Molecular Mass	71 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SDPEGEDWDE EAEDDGFDSD SSLSDSDLEQ DPEGLHLWNS FCSVDPYNPQ NFTATIQTAA RIVPEEPSDS EKDLGKSDL ENSSQSGSLP ETPEHSSGEE DDWESSADEA ESLKLWNSFC NSDDPYNPLN FKAPFQTSGE NEKGCGRDSKT PSESIVAISE CHTLLSCKVQ LLGSQESECP DSVQRDVLSG GRHTHVKRKK VTFLEEVTEY YISGDEDRKG PWEEFARDGC RFQKRIQETE DAIGYCLTFE

image



Background

Alternative Names Protein phosphatase 1 regulatory subunit 15B

Note

Product Name: Recombinant Human PPP1R15B Protein
Catalog #: PEH8930



For Research Use Only , Not for Diagnostic Use