


Product Name: Recombinant Human PPP1R42 Protein
Catalog #: PEH8931



Summary

Name	PPR42_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7Z4L9
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RSLKKLEKLY LGGNYIAVIE GLEGLGELRE LHVENQRLPL GEKLLFDPRT LHSLAKSLCI LNISNNNIDD ITDLELLENL NQLIAVDNQL LHVKDLEFLL NKLMKLVKID LNGNPVCLKP KYRDRILILVS KSLEFLDGKE IKNIERQFLM NWKASKDAKK ISKKRSSKNE DASNSLISNF KTMHHIVPVY YPQVGKPKLA FFSEIQRYPV

image



Background

Alternative Names	Protein phosphatase 1 regulatory subunit 42 (Leucine-rich repeat-containing protein 67)
--------------------------	---

Note

Product Name: Recombinant Human PPP1R42 Protein
Catalog #: PEH8931



For Research Use Only , Not for Diagnostic Use