


Product Name: Recombinant Human DNAAF6 Protein
Catalog #: PEH8932



Summary

Name	DAAF6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NQM4
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VDFESVSSVT ALEALSLLN PEEEDSDYG QTNGLSTIGA MGPGNIGPPQ IEELKVIPET SEENNEDIWN SEEIPEGAEY DDMWDVREIP EYEIIFRQQV GTEDIFLGLS KKDSSTGCCS ELVAKIKLPN

image



Background

Alternative Names Dynein assembly factor 6, axonemal (PIH1 domain-containing protein 3) (Sarcoma antigen NY-SAR-97)

Note

For Research Use Only , Not for Diagnostic Use