


Product Name: Recombinant Human PRRC2B Protein
Catalog #: PEH8934



Summary

Name	PRC2B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5JSZ5
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PREGTAFNIS SWDKNGSPNK QPSSEPEWTP EPRSSSSQHP EQTGRTRRS PIKKPVLKAL KVEDKEKELE KIKQELGEES TRLAKEKEQS PTAEKDEDEE NDASLANSSST TTLEDKGPGRH ATFGREATKF EEEEKPDKAW EARPPRESSD VPPMKRNNWI FIDEEQAFGV RGQARGRGRG FREFTFRGRP AGGNGSGLCG

image



Background

Alternative Names	Protein PRRC2B (HLA-B-associated transcript 2-like 1) (Proline-rich coiled-coil protein 2B)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PRRC2B Protein
Catalog #: PEH8934

