


Product Name: Recombinant Human SHROOM1 Protein
Catalog #: PEH8935



Summary

Name	SHRM1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q2M3G4
Host	E.coli
Species	Human
Molecular Mass	71 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GGPEPRTQSP GTDLLPYLDW DYVRVWGGP GPAPPDAALC TSPRPRPAVA ARSGPQPTEV PGTPGPLNRQ ATPLLYALAA EAEAAAQAAE PPSPPASRAA YRQRLQGAQR RVLRETSFQR KELRMSLPAR LRPTVPPARPP ATHPRSASLS HPGGEGEPAR SRAPAPGTAG RGPLANQQRK WCFSEPGKLD RVGRGGGPAR ECLGEACSSS GLPGPEPLEF QHPALAKFED HEVGWLPETQ PQGSMNLDSG SLKLGDAFRP ASRSRSASGE

image



Background

Alternative Names Protein Shroom1 (Apical protein 2)

Product Name: Recombinant Human SHROOM1 Protein
Catalog #: PEH8935



Note

For Research Use Only , Not for Diagnostic Use