


Product Name: Recombinant Human SMG9 Protein
Catalog #: PEH8937



Summary

Name	SMG9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H0W8
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LYGIERRRRW KEPGSGGPQN LSGPGGRERD YIAPWERERR DASEETSTSV MQKTPIILSK PPAERSKQPP PPTAPAAPPA PAPLEKPIVL MKPREEGKGP VAVTGASTPE GTAPPPPAAP APPKGEKEGQ RPTQPVYQIQ NRGMGTAAPA AMDPVVGQAK LLPPERMKHS IKLVDDQMNW CDSAIEYLLD QTDVLVVGVL GLQGTGKSMV MSLLSANTPE EDQRTYVFRA QSAEMKERGG

image



Background

Alternative Names Protein SMG9 (Protein SMG-9)

Note

Product Name: Recombinant Human SMG9 Protein
Catalog #: PEH8937



For Research Use Only , Not for Diagnostic Use