

Product Name: Recombinant Human SREK1IP1 Protein
Catalog #: PEH8938



Summary

Name	SR1IP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N9Q2
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VRAGCKKCGY PGHLTFECRN FLRVDPKRDI VLDVSSTSSE DSDEENEELN KLQALQEKRI NEEEEKKKEK SKEKIKLKKK

image



Background

Alternative Names	Protein SREK1IP1 (SFRS12-interacting protein 1) (SREK1-interacting protein 1) (Splicing regulatory protein of 18 kDa) (p18SRP)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use