


**Product Name: Recombinant Human TGM6 Protein**  
**Catalog #: PEH8941**



## Summary

<b>Name</b>	TGM3L_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	O95932
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	17 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	AHDTDQNLSV DKYVDSFGRT LEDLTEDSMW NFHVWNESWF ARQDLGPSYN GWQVLDATPQ EESEGVFRCG PASVTAIREG DVHLAHDGPF VFAEVNADYI TWLWHEDESR ERVYSNTKKI GRCISTKAVG

## image



## Background

**Alternative Names** Protein-glutamine gamma-glutamyltransferase 6 (EC 2.3.2.13) (Transglutaminase Y) (TGY) (TGase Y) (Transglutaminase-3-like) (TGase-3-like) (Transglutaminase-6) (TG6) (TGase-6)

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human TGM6 Protein**  
**Catalog #: PEH8941**

---

