


Product Name: Recombinant Human CHCHD2P9 Protein
Catalog #: PEH8945



Summary

Name	CHCH9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5T1J5
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RMAPPASRAP QMRAAPRPAP VAQPPAAAPP SAVGSSAAAP RQPGLMAQMA TTAAGVAVGS AVGHTQGHAV TGGFSGGSNA EPARPDIAVQ EPQGTQPAQQ

image



Background

Alternative Names	Putative coiled-coil-helix-coiled-coil-helix domain-containing protein CHCHD2P9, mitochondrial (Coiled-coil-helix-coiled-coil-helix domain-containing 2 pseudogene 9)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use