


Product Name: Recombinant Human MKRN4P Protein
Catalog #: PEH8946



Summary

Name	MKRN4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q13434
Host	E.coli
Species	Human
Molecular Mass	72 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QPYCGRTVPS CTEAPLQGSV TKEESEEQT AVETKKQLCP YAAVGQCRYG ENCVYLHGDL CDMCGLQVLH PMDAAQRSQH IQACIEAHEK DMEFSFAVQR SKDKVCGICM EVVYEKANPN EHRFGILSNC NHTFCLKCIR KWRSKEFES RIVKSCPQCR ITSNFVIPSE YWVEEKEEKQ KLIQKYKEAM SNKACKYFDE GRGSCPFGEN CFYKHMYPDG RREEPQRQQV GTSSRNPGQQ RNHFWEFFEE GANSNPFDDDE

image



Background

Alternative Names	Putative E3 ubiquitin-protein ligase makorin-4 (EC 2.3.2.27) (Makorin RING finger protein pseudogene 4) (Makorin RING finger protein pseudogene 5) (RING finger protein 64) (RING-type E3 ubiquitin transferase makorin-4) (Zinc
--------------------------	--

Product Name: Recombinant Human MKRN4P Protein
Catalog #: PEH8946



finger protein 127-Xp) (ZNF127-Xp) (Zinc finger protein 127-like 1)

Note

For Research Use Only , Not for Diagnostic Use