


Product Name: Recombinant Human GOLGA8F Protein
Catalog #: PEH8947



Summary

Name	GOG8F_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q08AF8
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ERVAGELQAQ VEYNQRISLL NEGQKERLRE QEERLQEQQE RLPEQEERLQ QLAEPQNSFK ELNNENKSVL QLEQQVKELQ EKLGKERLEA ASQQKQQLTA QLSLMALPGE GDGGGHLDSE GEEAPRPIPS IPQDLESREA MSGFMDHLEE KADLSELVEK EELGFFQYYR ERCHQKVYHP ITKPGGSAKD AAPGGGHHQA GPGQGGDEGE AAGAAGDGVA AGGDYKGHSK FLVTAQNPAH EPSPGAPAPQ ELGAAHKHGD LCEVSLTDSV EPVQGEAREG SPHDNPTAQP IVQDHQEHGP

image



Background

Alternative Names Putative golgin subfamily A member 8F/8G

Product Name: Recombinant Human GOLGA8F Protein
Catalog #: PEH8947



Note

For Research Use Only , Not for Diagnostic Use