


Product Name: Recombinant Human HSP90AA4P Protein
Catalog #: PEH8948



Summary

Name	HS904_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q58FG1
Host	E.coli
Species	Human
Molecular Mass	49 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MESLTDPSKL DSGKEPHISL IPNKQDRTLTV IVDTGIGMTK ADLINNLGTI TKSETKVFME VLQAGADISM IGQFSVGFYS AYSVAEKVTV ITKHNNDEQY AWESSLRGSF TEYREFYKSL TINWEDYLAV KHFSVEGQLE FRAFLFVPRL APFELLETRK KKNKIKLSAR RDLIMDNCEE LIPEYLNfir GVVDSIDLPL NIFRETKDQV ANSTIVQRLW KHGLEVIYTI EPIDEYCVQQ LKEFEGKTLV SVTKEDLELP EDEEEKKKQE EGKQKTKQKK NQSLRTSAKS TYGWTANMER IMKAQALRDN STTGYYMAAKK HLEINPDHSF IDTLRQKAET DKNDKSVKDL VILLYETALL SSDFGLEGPQ THANRIYRMN KLGLGTDEDD PTADDTSAAV TEEMPPLEGD DDTSRMEK

image



Background



Product Name: Recombinant Human HSP90AA4P Protein
Catalog #: PEH8948

Alternative Names

Putative heat shock protein HSP 90-alpha A4 (Heat shock 90 kDa protein 1 alpha-like 2) (Heat shock protein 90-alpha D) (Heat shock protein 90Ad)

Note

For Research Use Only , Not for Diagnostic Use