


Product Name: Recombinant Human HSP90AB2P Protein
Catalog #: PEH8949



Summary

Name	H90B2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q58FF8
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NLGTIAKFQD QTEYLEEMQV KEVVEKHSQF LGYPITLYLE KEREKEISDG KAEEEKGEKE EENKDDEEKP KIEDVGSDEE DDSGKD KKKK TKKIKEYID QEELNKT KPI WTRNTEDITQ EEYGEFYKSL TNDWKDHLAV RYFSVEEYVS RMKEIQKSIY YITGESKEQV ANSAFVEQVW KRDSRVVYMT EPIDGYQLKE FDGKSLVSVT KEGLELPEDG EEKKRMEERK AKFENLCKFM KETLDKKVEM

image



Background

Alternative Names	Putative heat shock protein HSP 90-beta 2 (Heat shock protein 90-beta b) (Heat shock protein 90Bb)
--------------------------	--

Note



Product Name: Recombinant Human HSP90AB2P Protein
Catalog #: PEH8949

For Research Use Only , Not for Diagnostic Use