


Product Name: Recombinant Human NSUN7 Protein
Catalog #: PEH8951



Summary

Name	NSUN7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NE18
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FLSILTRERD PSETVSVNDV LARAAAKGLL DGIELGKSSK REKKKKKSKT SLTKGATTDN GIQMKIAEFL NRETKASANL SETVTKPPLP QKNTAQVGAS SQTRKPNKLA PHPAVPAFVK

image



Background

Alternative Names	Putative methyltransferase NSUN7 (EC 2.1.1.-) (NOL1/NOP2/Sun domain family member 7)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use