


Product Name: Recombinant Human PMS2P1 Protein
Catalog #: PEH8952



Summary

Name	PM2P1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A4D2B8
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KERREVFPP RPERVEHNVE SSRWEPRRRG ACGSRGGNFP SPRGGSGVAS LERAESSSTE PAKAIKPIDR

image



Background

Alternative Names	Putative postmeiotic segregation increased 2-like protein 1 (PMS2-related protein 2) (Postmeiotic segregation increased 2-like protein 13) (Postmeiotic segregation increased 2-like protein 6) (Postmeiotic segregation increased 2-like protein 8) (Postmeiotic segregation increased protein 3) (hPMS3) (Postmeiotic segregation increased protein 8) (Putative postmeiotic segregation increased 2 pseudogene 1)
--------------------------	--

Product Name: Recombinant Human PMS2P1 Protein
Catalog #: PEH8952



Note

For Research Use Only , Not for Diagnostic Use