

Product Name: Recombinant Human ST13P4 Protein
Catalog #: PEH8953



Summary

Name	ST134_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IZP2
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AFVKMCKKDP SILHTQEMRF LREWVESMGG TATQKAKSEE NTKEEKPDSK VEEDLKADEP SSEESDLEID KEGVIEPDTD APQEMGDENA EITEEVMDQA NDKKVAAIEA LNDGELQKAI DLFTDAIKLN PRLAILYAKR ASVFVKLQKP NAAIRDCDRA IEINPDSAQP YKRRGKAHRL LGHWEEAAHD

image



Background

Alternative Names Putative protein FAM10A4 (Suppression of tumorigenicity 13 pseudogene 4)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ST13P4 Protein
Catalog #: PEH8953

