


Product Name: Recombinant Human RABL2B Protein
Catalog #: PEH8956



Summary

Name	RBL2B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UNT1
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAEDKTKPSE LDQGKYDADD NVKIICLGDS AVGKSKLMER FLMDGFQPPQ LSTYALTLYK HTATVDGRTI LVDFWDTAGQ ERFQSMHASY YHKAHACIMV FDVQRKVTYR NLSTWYTELR EFRPEIPCIV VANKIDDINV TQKSFNFAKK FSLPLYFVSA ADGTNVVKLF NDAIRLAVSY KQNSQDFMDE IFQELENFSL EQEEEDVPDQ EQSSSIETPS EEAASPHS

image



SDS

Background

Alternative Names Rab-like protein 2B

Note

Product Name: Recombinant Human RABL2B Protein
Catalog #: PEH8956



For Research Use Only , Not for Diagnostic Use