


Product Name: Recombinant Human RSPH6A Protein
Catalog #: PEH8957



Summary

Name	RSH6A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H0K4
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	INCDLSLYEH LVNLLTKILN QRPEDPLSVL ESLNRTTQWE WFHPKLDTLR DDPEMQPTYK MAEKQKALFT RSGGGTEGEQ EMEEEVGETP VPNIMETAFY FEQAGVGLSS DESFRIFLAM KQLVEQQPIH TCRFWGKILG IKRSLVAEV EFREGEEAE EEEVEEMTEG GEVMEAHGEE EGEEDEEKAV DIVPKSVWKP PPVIPKEESR SGANKLYFV CNEPGLPWTR LPHVTPAQIV NARKIKKFFT GYLDTPVVS Y PFPGNEANY LRAQIARISA ATQVSPPLGFY QFSEEEGDDEE EEGGAGRDSY EENPDFEGIP VLELVDSMAN

image



Background

Alternative Names	Radial spoke head protein 6 homolog A (Radial spoke head-like protein 1)
--------------------------	--

Product Name: Recombinant Human RSPH6A Protein
Catalog #: PEH8957



Note

For Research Use Only , Not for Diagnostic Use