


Product Name: Recombinant Human REPS2 Protein
Catalog #: PEH8960



Summary

Name	REPS2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His
Accession #	Q8NFH8
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PPTLQPEYLQ AAFPKPKWDC QLFDSEYSEL PANQQPRDLN RMEKTSVKDM ADLPVPNQDV TSDDKQALKS TINEALPKDV SEDPATPKDS NSLKARPRSR SYSSTSIEEA MKRGEDPPTP PPRPQKTHSR ASSLDLNKVF QPSVPATKSG LLPPPPALPP RPCPSQSEQV SEAELLPQLS RAPSQAAESE PAKKDVLYSQ PPSKPIRRKF RPENQATENQ EPSTAASGPA SAATMKPHPT VQKQSSKQKK AIQTAIRKNK EANAVALRLN SELQQQLKEV HQERIALENQ LEQLRPVTVL

image



Background

Alternative Names	RalBP1-associated Eps domain-containing protein 2 (Partner of RalBP1) (RalBP1-interacting protein 2)
--------------------------	---

Product Name: Recombinant Human REPS2 Protein
Catalog #: PEH8960



Note

For Research Use Only , Not for Diagnostic Use