


Product Name: Recombinant Human RGP4 Protein
Catalog #: PEH8961



Summary

Name	RGPD4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7Z3J3
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KEGFSIPVSA DGFKFGISEP GNQEKESKP LENDTGFAQQ DISGQKNGRG VIFGQTSSTF TFADVAKSTS GEGFQFGKKD PNFKGFSGAG EKLFSQCGK MANKANTSGD FEKDDDAYKT EDSDDIHFEF VVQMPEKVEL VIGEEGEKVL YSQGVKLFRR DAEVRQWKER GLGNLKILKN EVNGKPRMLM RREQVLKVCA NHWITTTMNL KPLSGSDRAW MWSASDFS DG DAKLERLAAK FKTPELAEFF KQKFECCQQL LLDIPLQTPH KLVDTGRAAK LIQRAEEMKS GLKDFKTFLT NDQTKVTEEE NKGSGTGAAG

image



Background

Alternative Names RanBP2-like and GRIP domain-containing protein 4

Product Name: Recombinant Human RGP4 Protein
Catalog #: PEH8961



Note

For Research Use Only , Not for Diagnostic Use