


## Summary

<b>Name</b>	RGPD5_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q99666
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	38 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	SSNTEFKSTK EGF SIPVSAD GFKFGISEPG NQEKKREKPL ENDTGFQAQD ISGRKKGRGV IFGQTSSTFT FADVAKSTSG EGFQFGKKDL NFKGFSGAGE KLFSSRYGKM ANKANTSGDF EKDDDAYKTE DSDDIHFEPV VQMPEKVELV TGEEGEKVLV SQGVKLFKRFD AEVRQWKERG LGNLKILKNE VNGKLRMLMR REQVLKVCAN HWITTTMNLK PLSGSDRAWM WSASDFSDGD AKLERLAAKF KTPELAEFEK QKFEECQRL LDIPLQTPHK LVDTGRAAKL IQRAEEMKSG LKDFKTLTN DQTKVTEEN KGSGTGAAGA

## image



## Background

<b>Alternative Names</b>	RANBP2-like and GRIP domain-containing protein 5/6 (Ran-binding protein 2-like 1/2) (RanBP2-like 1/2) (RanBP2L1) (RanBP2L2) (Sperm membrane
--------------------------	---

**Product Name: Recombinant Human RGD5 Protein**  
**Catalog #: PEH8962**



---

protein BS-63)

**Note**

For Research Use Only , Not for Diagnostic Use