


Product Name: Recombinant Human RCC1L Protein
Catalog #: PEH8969



Summary

Name	RCC1L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96I51
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SAACGYGFTL LSKTADVTK VWGMGLNKDS QLGFRSRKD KTRGYEYVLE PSPVSLPLDR PQETRVLQVS CGRAHSLVLT DREGVFSMGN NSYGQCGRKV VENEIYSESH RVHRMQDFDG QVVQVACGQD HSLFLTDKGE VYSCGWGADG QTGLGHYNIT SSPTKLGGDL AGVNVIQVAT YGDCCLAUSA DGGLFGWGNS EYLQLASVTD STQVNVPRCL HFSGVGKVRQ AACGGTGCAV LNGEHVFVW GYGILGKGPV LVESAVPEMI PPTLFGLTEF NPEIQVSRIR CGLSHFAALT NKGELFVWGK NIRGCLGIGR LEDQYFPWRV TMPGEPVDVA CGVDHMTLA

image



Background

Alternative Names	RCC1-like G exchanging factor-like protein (RCC1-like) (Williams-Beuren syndrome chromosomal region 16 protein)
--------------------------	---

Product Name: Recombinant Human RCC1L Protein
Catalog #: PEH8969



Note

For Research Use Only , Not for Diagnostic Use