


Product Name: Recombinant Human RHOV Protein
Catalog #: PEH8977



Summary

Name	RHOV_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q96L33
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MPPRELSEAE PPPLRAPTPP PRRRSAPPEL GIKCVLVGDG AVGKSSLIVS YTCNGYPARY RPTALDTFSV QVLVDGAPVR IELWDTAGQE DFDRLRSLCY PDTDVFLACF SVVQPSSFQN ITEKWLPEIR THNPQAPVLL VGTQADLRDD VNVLIQLDQG GREGVPVQPQ AQGLAEKIRA CCYLECSALT QKNLKEVFDS AILS AIEHKA RLEKKLNAKG VRTL SRCRWK KFFCFV

image



Background

Alternative Names	Rho-related GTP-binding protein RhoV (CDC42-like GTPase 2) (GTP-binding protein-like 2) (Rho GTPase-like protein ARHV) (Wnt-1 responsive Cdc42 homolog 2) (WRCH-2)
--------------------------	--

Product Name: Recombinant Human RHOV Protein
Catalog #: PEH8977



Note

For Research Use Only , Not for Diagnostic Use