


Product Name: Recombinant Human RPS6KC1 Protein
Catalog #: PEH8979



Summary

Name	KS6C1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96S38
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FRIDSKDSAS ELLGLDFGEK LYSLKSEPLK PFFTLPDGDS ASRSFNTSES KVEFKAQDTI SRGSDDSVPV ISFKDAAFDD VSGTDEGRPD LLVNLPGELE STREAAAMGP TKFTQTNIGI IENKLLEAPD VLCLRLSTEQ CQAHEEKGIE ELSDPSGPKS YSITEKHYYAQ EDPRMLFVAA

image



Background

Alternative Names	Ribosomal protein S6 kinase delta-1 (S6K-delta-1) (EC 2.7.11.1) (52 kDa ribosomal protein S6 kinase) (Ribosomal S6 kinase-like protein with two PSK domains 118 kDa protein) (SPHK1-binding protein)
--------------------------	--

Note

Product Name: Recombinant Human RPS6KC1 Protein
Catalog #: PEH8979



For Research Use Only , Not for Diagnostic Use