


Product Name: Recombinant Human RILPL2 Protein
Catalog #: PEH8981



Summary

Name	RIPL2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q969X0
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DEVGPEGALG KSPFQLTAED VYDISYLLGR ELMALGSDPR VTQLQFKVVR VLEMLEALVN EGSLALEELK MERDHLRKEV EGLRRQSPPA SGEVNLGPNK MVVDLTDPNR PRFTLQELRD VLQERNKLKS QLLVVQEELQ CYKSGLIIPPR EGPGGRREKD AVVTSAKNAG RNKEEKTIK

image



Background

Alternative Names	RILP-like protein 2 (Rab-interacting lysosomal protein-like 2) (p40phox-binding protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RILPL2 Protein
Catalog #: PEH8981

