


Product Name: Recombinant Human RNF113B Protein
Catalog #: PEH8982



Summary

Name	R113B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IZP6
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DVVYRSTRSA KPVGPEDMGA TADFEQDTEK EHHTPTILKC SQRVQEALRG REHDHIYRGI HSYLRYLKP K DTSMGNSSSG MARKGPIRAP GHLRATVRWD YQPDICKDYK ETGFCGFGDS CKFLHDRSDY KLGWEIEREL EEGRYCICED ENHEVGSEEE EIPFRFCFICR QAFQNPVVTK

image



Background

Alternative Names RING finger protein 113B (Zinc finger protein 183-like 1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RNF113B Protein
Catalog #: PEH8982

