



Product Name: Recombinant Human RNF207 Protein
Catalog #: PEH8983



Summary

Name	RN207_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZRF8
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SLQLDVQIAS EHASLEGMRV VFQEIWEEAY QRVANEQEIY EAQLHDLLQL RQENAYLTTI TKQITPYVRS IAKVKERLEP RFQAPVDEQS ESLQNTHTDD RNNAASARNN PGSVPEKREK TSEPKGNSWA PNGLSEEPLL

image



Background

Alternative Names RING finger protein 207

Note

For Research Use Only , Not for Diagnostic Use