


Product Name: Recombinant Human RADX Protein
Catalog #: PEH8987



Summary

Name	RADX_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6NSI4
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DLQYREQKRI AIQGIITAIK YIPHSSATES ASASETLRNA NRPSTSQAAR VEIQERNGKR HQDDEPVNSQ YFQTTSTNLS LSNKIRILQG PHANPVAVPQ PGASVQTKGI KPGMPSIFNR

image



Background

Alternative Names RPA-related protein RADX (RPA-related and RAD51-antagonist, X-chromosome)

Note

For Research Use Only , Not for Diagnostic Use