



Product Name: Recombinant Human MRM3 Protein
Catalog #: PEH8988



Summary

Name	MRM3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His
Accession #	Q9HC36
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VFPSGEVVEQ KRAPGKQPRK APSEASAEQ REKQPLEESA SRAPSTWEES GLRYDKAYPG DRRLLSSVMTI VKSRPFREKQ GKILLEGRRL ISDALKAGAV PKMFFFSRLE YLKELPVDKL KGVSLIKVKF EDIKDWSDLV TPQGIMGIFA KPDHVKMTYP KTQLQHSLPL LLICDNL RDP GNLGTILRSA AGAGCSKVLL TKGCVDAWEP KVL RAGMGAH FRMPIINNLE WETVPNYLPP DTRVYVADNC GLYAQAEMSN KASDHGWVCD QRVMKFH KYE EEEDVETGAS QDWLPHVEVQ SYDSDWTEAP AAVVIGGETY GVSLESLQLA ESTGGKRLLI PVVPGVDSL N SAMAASILLF EGKRQLRGRA EDLSRDRSYH

image



Background

Alternative Names rRNA methyltransferase 3, mitochondrial (EC 2.1.1.-) (16S rRNA

Product Name: Recombinant Human MRM3 Protein
Catalog #: PEH8988



(guanosine(1370)-2'-O)-methyltransferase) (16S rRNA [Gm1370] 2'-O-methyltransferase) (RNA methyltransferase-like protein 1)

Note

For Research Use Only , Not for Diagnostic Use