


Product Name: Recombinant Human SRBD1 Protein
Catalog #: PEH8989



Summary

Name	SRBD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N5C6
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QVQDVVLKDE FSSFSELSSA SEEDDKEDSA WEPQKKVPRS RKQPPPKESK PKRMPRVKKN APQISDGSEV VVVKEELNSS VAIADTALED RKNKLDTVQT LKTAKTKQKC AAQPHTVRRT KKLKVEEETS KASNLEGESN SSETPSTSTV WGGTCKKEEN DDDFTFGQSA LKKIKTETYP QGQPVKFPAN ANSTKEEVEM NWDMMVQVLSE RTNIEPWVCA NIIRLFNDDN TIPFIIRYRK ELINNLDADS LREVQQTLEE LRVAKKVHS TIQKIKKEGK MSECLLKAML NCKTFEELEH VSAPYKTGSK GTKAQRARQL GLEGAARALL EKPGEISLLS YIRPDVKGLS TLQDIEIGVQ HILADMIKD KDTLDFIRNL

image



SDS

Background

Alternative Names	S1 RNA-binding domain-containing protein 1
--------------------------	--

Product Name: Recombinant Human SRBD1 Protein
Catalog #: PEH8989



Note

For Research Use Only , Not for Diagnostic Use